


Brain tissue preparation and immunohistochemistry

RS Renato Socodato TS Teresa Summavielle JR João B. Relvas

Updated date: Oct 23, 2020

 An abbreviated version of this protocol was published in Science Signaling in Sep 2020

Daily alcohol intake triggers aberrant synaptic pruning leading to synapse loss and anxiety-like behavior

DOI: 10.1126/scisignal.aba5754

Detailed protocol

Brain tissue preparation and immunofluorescence for evaluating synaptic pruning by microglia

Joana Tedim-Moreira^{1,*}, Tiago O. Almeida^{1,*}, Camila C. Portugal¹, Teresa Summavielle¹, João B. Relvas^{1,2} and Renato Socodato¹

¹Instituto de Investigação e Inovação em Saúde and Instituto de Biologia Molecular e Celular, Universidade do Porto, 4200-135 Porto, Portugal²Department of Biomedicine, Faculty of Medicine, University of Porto, 4200-319 Porto, Portugal.

*These authors contributed equally to this protocol.

Email to: renato.socodato.ibmc.up.pt or renato.socodato@gmail.com

Step-by-step:

• Tissue preparation:

- 1) Perfuse the mouse transcardially with ice-cold PBS (~10ml) and then with 4% PFA (~25-30ml) using an osmotic pump.
 - 2) Dissect the brain tissue (with the olfactory bulb and cerebellum) and postfix it in 4% PFA in PBS, pH 7.2, for 12-18h.
 - 3) Wash the tissue once with excess PBS.
 - 4) Transfer the tissue to a 15 ml tube and cryoprotect the tissue using a sucrose gradient in a row — 15% sucrose then 30% sucrose (24h in each solution or until the tissue reaches the bottom of the tube).
 - 5) After 24h, embed the tissue in optimal cutting temperature embedding medium (OCT – [cat# 6769006; ThermoFisher](#)) and store at -20°C (short-term) or -80°C (long-term).
 - 6) Using a CM3050S Cryostat (Leica Biosystems) (chamber temperature at -24°C and the object temperature at -18°C), cut 30µm coronal sections of the desired brain regions.
 - 7) Carefully collect the brain slices, non-sequentially, on Superfrost Plus slides ([cat# J1800AMNZ; ThermoFisher](#)).
- Tissue sections encompassing identical stereological regions from the brains of different mice should be collected on the same glass slides and stored at -20°C.
- 8) Thaw slides with desired tissue sections for at least 1h before starting the immunostaining procedure.

• Immunostaining:

- 1) Hydrate the sections with PBS for 15 min.
- 2) Permeabilize the sections for 15 min.
- 3) Remove the permeabilization solution and wash sections in excess PBS for 10 min.
- 4) Carefully remove excess PBS and add the blocking solution for 1h at 22-25°C.
- 5) Remove the blocking solution and incubate with primary antibodies (prepared in blocking solution) in a humidified chamber at 4°C for 12-18h.
- 6) Remove primary antibodies and wash sections in excess PBS 3x10 min.
- 7) Remove PBS and incubate sections with secondary antibodies (1:600 AlexaFluor 488 [[cat# A110081](#)], 1:400 AlexaFluor 568 [[cat# A110771](#)], and 1:500 AlexaFluor

647 [cat# A21325]; all from ThermoFisher) in blocking solution for 2h at 22-25°C.

➤ Depending on the microscope system available, antibodies` titration might require further optimizations.

8) Remove secondary antibodies and wash sections in excess PBS 3x10 min.

➤ If performing double or triple-labeling immunostaining, add the primary antibody followed by the adequate secondary antibody and repeat the process for each combination of primary/secondary antibodies. Never combine the primary antibodies in the same incubation step with a single master mix.

9) Coverslip the slides using Fluoroshield anti-fading reagent (cat# F6182; Sigma) and store at 4°C.

10) Visualize the slides under a Leica TCS SP8 confocal microscope.

Solutions:

1. Permeabilization solution: 0.25% TritonX-100 in PBS

2. Blocking solution: 5% BSA (cat# MB04602; Nzytech) + 5% FBS (cat# 10270106; ThermoFisher) + 0.1% TritonX-100 in PBS

Tested primary antibodies:

Antibody	Dilution	Company	Catalogue Number
Iba1	1:500	Wako	019-19741;RRID:AB_839504
vGlut-1	1:1,000	Synaptic Systems	135 303; RRID:AB_887875
PSD95 -clone 6G6-1C9	1:600	Thermo Scientific	catalog no. MA1-045;RRID:AB_325399
CD68 - clone FA-11	1:400	Bio-Rad	MCA1957T; RRID:AB_2074849

How to cite:(Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Socodato, R. , Summavielle, T. and Relvas, J. (2020). Brain tissue preparation and immunohistochemistry. Bio-protocol Preprint. bio-protocol.org/prep567.
2. Socodato, R., Henriques, J. F., Portugal, C. C., Almeida, T. O., Tedim-Moreira, J., Alves, R. L., Canedo, T., Silva, C., Magalhães, A., Summavielle, T. and Relvas, J. B.(2020). Daily alcohol intake triggers aberrant synaptic pruning leading to synapse loss and anxiety-like behavior . Science Signaling 13(650). DOI: [10.1126/scisignal.aba5754](https://doi.org/10.1126/scisignal.aba5754)

Copyright: Content may be subjected to copyright.